




3 KEYS TO SUCCESS WITH EXOS

& caring for your skin and brace


1. Not Too Tight – Wear it Right and leave some “WIGGLE ROOM”

- Exos Braces are meant to be worn under light pressure and not tight to the skin. A little wiggle room inside an Exos Brace promotes air circulation and helps maintain dry healthy skin.
- An over-tightened brace creates “shear” or pressure on the skin during motion, and can result in irritation, rash, odor, and skin issues.



2. If the brace and skin get wet, they “CAN’T STAY WET”


- The Brace may be washed or worn while swimming or bathing, only if the physician allows patient to loosen or remove the brace. If not, the patient cannot get the brace and skin wet.
- It is imperative that the inside surface of the brace and the patient’s skin underneath are completely dry after swimming or bathing.
- See below for cleaning and drying instructions.



3. Keep your “SKIN AND BRACE CLEAN AND DRY”

- If your physician has given permission to loosen or remove the brace, wash the brace (inside and out) and your skin with anti-microbial soap, rinse thoroughly, and completely dry the skin and brace (inside and out) with a hair dryer.
- A clean and dry environment inside the brace will help maintain healthy skin and reduce odor and skin issues.

1. Loosen or remove the brace as instructed by physician.
2. Wash inside and outside of brace and skin with antibacterial soap and water.
3. Rinse thoroughly with water.
4. Use a hair dryer set on high volume and cool setting to thoroughly dry the interior of the brace and skin.
5. Retighten the brace.



Warnings:

- Consult your physician if unusual swelling, increasing discomfort, pain, numbness or discoloration appears.
- Instruction for wearing, opening, closing or removing the brace should be advised by a physician.
- If brace does not fit properly, or is irritating, it may be reformed and adjusted by a health professional.
- If brace breaks or cracks, see your health professional for a replacement.
- Brace should be adjusted so that is supportive, but not constrictive.
- There should be a slight amount of space under the brace for air to circulate, and keep the skin dry and healthy.
- Wearing the brace too tight may result in skin irritation, odor, or other skin issues.
- Wearing the brace too loose will reduce support, and may compromise healing.
- Completely dry skin underneath brace after getting wet.
- Do not expose brace to moisturizers, chemicals, or solvents, as they may affect durability or cause skin issues.
- Do not expose brace to heat over 130°F (55°C)



DJO, LLC.
1430 Decision Street
Vista, CA 92081-8553 U.S.A.
Phone 760.727.1280
Toll Free 800.526.8785
Fax 800.457.4221